







Compact Mode



Trolley Mode



Split Mode



The best travel mobility scooter

DIMENSIONS	WEIGHT	
Trolley mode H 28" / 71 cm L 15.3" / 38.9 cm W 16.5" / 41.9 cm	Front half 26.4 lbs / 12 kg	Rear half 37.5 lbs / 17 kg
Drive mode H 35" / 89 cm L 47" / 120 cm W 22" / 56 cm	Assembled 64 lbs / 29 kg (w/o battery)	Battery 4.1 lbs / 1.85 kg
WEIGHT CAPACITY	SPEED	
Maximum user weight 264 lbs / 120 kg	Most countries, 4 mph / 6.4 km per hour*	

BATTERY

48 Volts Lithium-ion, more than 20 km / 12 miles of continuous driving, only 4 hours on average for full charge, charge on-board and off-board

COMPATIBILITY



Airline compatible



Travel easy



Great for cruises









MORE FOR ZITTLE

ATTO ACCESSORIES are created specifically to enhance, support and work with the ATTO mobility scooter in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.




















ATTO Left Hand Throttle



ATTO Standard Flight Battery



ATTO XL Flight Battery

